



ratiolab® Accupetta Electronic Pipetting Aid

CE IVD

- can fill 25 ml pipet in less than 4 seconds
- safety valve and hydrophobic filters provide double protection against fluid penetration
- two different speed modes – high & low along with gravity drain
- the pump speed can also be fine-tuned by varying finger pressure on operating knobs for better control of speed
- cadmium-free environment-friendly NiMH batteries
- intelligent charger prevents over-charging/heating of batteries
- low battery indicator, power 2 x 1.2 V/1000 mAh
- battery service life: 8 hours continuous use
- battery charging time: 2 – 3 hours
- volumetric range: plastic or glass pipets from 0.1 ml to 100 ml
- autoclavable pipet cone I
- 1 year guarantee

Order No. 32 00 300 Accupetta contains:

- 1x Accupetta
- 2x Rechargeable battery
- 1x Adapter for battery charger
- 2x Spare filter 0,45 µm
- 1x Wall mount for Accupetta

Suitable serological pipets see on page 40.

| Product | Quantity per Pack | Order No. |
|--|-------------------|-----------|
| ratiolab® Accupetta | 1 | 32 00 300 |
| Accessories: | | |
| Spare silicone pipet adapter | 1 | 32 00 301 |
| Spare filter hydrophobic, non-sterile, 0.2 µm | 1 | 32 00 302 |
| Spare filter hydrophobic, non-sterile, 0.45 µm | 1 | 32 00 303 |
| Battery for Accupetta, set of 2 | 1 | 32 00 306 |
| Battery charger for Accupetta | 1 | 32 00 309 |
| Adapter for battery charger | 1 | 32 00 311 |
| Wall mount for Accupetta | 1 | 32 00 312 |