

# Multi Purpose Beakers and Multi Purpose Beaker System

## Multi Purpose beakers with screw-caps

CE IVD

Polypropylene (PP) or polystyrene (PS).

For sampling, transport and storage of solid, pasty and liquid samples.

- with red screw-cap
- supplied in non-sterile or sterile, individually packed



Volume ml	Fig. No.	Material	Graduation and Labeling area	Sterilized	Height with cap, mm	Ø cap, outer, mm	Packaging	Quantity per Pack	Order No.
30	1	PP	–	–	47	43	cardboard box 1 x 1000	1000	59 21 173
30	1	PP	–	■	47	43	bags 500 x 1	500	59 21 174
60	2	PP	■	–	66	43	cardboard box 1 x 600	600	59 21 175
60	2	PP	■	■	66	43	bags 500 x 1	500	59 21 176
120	3	PP	■	–	73	63	cardboard box 1 x 500	500	59 21 182
120	3	PP	■	■	73	63	bags 250 x 1	250	59 21 183
200	4	PP	–	–	80	64	cardboard box 1 x 200	200	59 21 185
200	4	PP	–	■	80	64	bags 150 x 1	150	59 21 184
200	5	PS	–	■	82	64	bags 150 x 1	150	59 21 188

## Multi Purpose Beakers

Glass-clear polystyrene (PS).

- volume 100 ml, with graduation
- self-adhesive film cover, leakproof and easily labelled
- comfortable sampling, the film cover can be pierced with pipet tips
- practical packaging, 50 multi purpose beakers are stacked in PE-bags for use with the beaker dispenser



Product	Packaging	Quantity per Pack	Order No.
ratiolab® multi-purpose beakers with cover	bags 20 x 50	1000	59 21 170
ratiolab® multi-purpose beakers without cover	bags 20 x 50	1000	59 21 171
self-adhesive film covers for beakers	strip 125 x 4	500	59 21 160

## Beaker Tray and Beaker Dispenser for Multi Purpose Beakers

Beaker Tray: Robust plastics with easy-to-clean surface.

- for 10 ratiolab® multi purpose beakers

Beaker Dispenser: Robust plastics with easy-to-clean surface.

- for 50 ratiolab® multi purpose beakers
- for individual beaker dispensing
- glass-clear dispenser door for refills makes stock visible
- separate side box for stocking of self-adhesive film covers



Product	Fig. No.	Quantity per Pack	Order No.
ratiolab® beaker tray	1	1	59 21 200
ratiolab® beaker dispenser	2	1	59 21 190

# Multi Purpose bottles, Petri Dishes and Inoculating Loops



## Multi Purpose bottles

Polyethylene, blue translucent (PE).

- inclusive lid with nonwoven fabric liner
- with graduations, light protection
- wide neck diameter 80 mm
- Volume: 2.5 l

Product	Packaging	Quantity per Pack	Order No.
Multi purpose bottle 2.5 l with self-adhesive labels	cardboard box 1 x 60	60	59 22 250
lids for Multi purpose bottle 2.5 l	bag 1 x 100	100	59 22 204
self-adhesive labels for Multi purpose bottle 2.5 l	bag 1 x 100	100	59 22 253



## Petri Dishes

CE IVD

Transparent polystyrene (PS).

- uncoated, with vents
- sterilized or clean as produced

Ø* x H mm	Description	Packaging	Quantity per Pack	Order No.
60 x 14	with vents	50 x 10	500	90 50 006
60 x 14	sterilized, with vents	50 x 10	500	90 50 106
90 x 15	with vents	25 x 20	500	90 50 000
90 x 15	sterilized, with vents	25 x 20	500	90 50 100
147 x 18	with vents	11 x 20	220	90 50 004
147 x 18	sterilized, with vents	11 x 20	220	90 50 104
65 x 14	with triple vents with compartments	25 x 20	500	90 50 005
65 x 14	sterilized, with triple vents with compartments	25 x 20	500	90 50 105

\*external diameter



## Inoculating Loops

CE IVD

Made of polystyrene (PS).

- hexagonal handle for safe handling
- beta-sterile
- sterile, individually packed in easy to open peel-packs

Product	Packaging	Quantity per Pack	Order No.
1 µl inoculating loops	100 peel-packs of 20	2000	70 20 100
10 µl inoculating loops	100 peel-packs of 20	2000	70 20 200