



## Cover Films

For micro test plates, deep-well plates, micro test tubes and PCR plates  
Adhesive base must be free of grease, dust, degassing and air-conditioned in order to ensure safe bonding.

- self-adhesive
- for manual operation
- adhesive with high solvent resistance
- all films of SBS standard format

### 1 Soft-PVC peeling cover film

- self-adhesive with one or two tabs
- low-cost film for polystyrene plates
- can be removed without any residues
- minimum application temperature: + 5 °C
- temperature resistant from – 20 °C to + 70 °C
- ideal for the storage of plates at room temperature

### 2 Aluminum universal cover film

- can be easily removed
- ideal for piercing
- for all plate types; PS, PP, PC, PE
- temperature resistant from – 80 °C to + 120 °C
- ideal for the storage of plates at room temperature (20 °C-23 °C) and dry room

### 3 Piercing-foil with cavity circle 96-well

- each cavity is marked by a circle and a numerical printing
- prepierced in each cavity optimizes opening by tip
- self-closing film reduces the risk of cross-contamination
- also for use with micro-racks and critical PCR plates
- temperature resistant from – 80 °C to + 30 °C
- DMSO-resistant

### 4 Gas permeable peeling film

- SBS-Standard
- liquid tight
- high permeability
- temperature resistant from – 5 °C to + 60 °C
- ideal for the storage of plates at room temperature (20 °C-23 °C) and dry room

### Polypropylene peeling cover film

- suitable for Real-Time-PCR
- self-adhesive with one or two tabs
- can be removed without any residues and repeatedly used
- for all plate types; PS, PP, PC, PE, DMSO-resistant
- temperature resistant from – 80 °C to + 120 °C
- ideal for the storage of plates at room temperature (20 °C-23 °C) and dry room

Further films on request.

Product	Fig No.	Material	Packaging	Quantity per Pack	Order No.
soft-pvc peeling cover film	1	PVC	pack 1 x 100	100	60 18 610
aluminum universal cover film	2	AL	pack 1 x 100	100	60 18 611
piercing-foil with cavity circle 96-well	3	PP	pack 1 x 100	100	60 18 619
gas permeable peeling film	4	PP	pack 1 x 50	50	60 18 618
polypropylene peeling cover film	-	PP	pack 1 x 100	100	60 18 612
polyvinyl chloride peeling cover film with 2 fixing bracket/tearing slices	-	PVC	pack 1 x 100	100	60 18 614
polypropylen peeling cover film with 2 fixing bracket/tearing slices	-	PP	pack 1 x 100	100	60 18 620

PVC = Polyvinyl chloride, AL = Aluminium, PP = Polypropylene