

# Micro Test Plates

## 96-Well Micro Test Plates, volume approx. 300 µl

Clear polystyrene (PS) or translucent polypropylene (PP).

- with U-, F-, or V-bottom
- SBS standard micro test plate format
- 8 x 12 wells alphanumerically coded
- wells with raised rims for reduced cross-contamination and effective sealing with Ratiolab cover films (see on page 69)
- precise design, perfectly shaped for use with automatic systems
- optimum space for precise positioning
- space-saving and stackable
- polypropylene plates may be sealed with heat-sealing films without deformation



### Volume approx. 300 µl

Bottom	Fig. No.	Material	Sterilized	Packaging	Quantity per Pack	Order No.
U	1	PS	–	bags 10 x 10	100	60 18 111
U	1	PS	■	individually in bags 100 x 1	100	60 18 113
U	2	PP	–	bags 10 x 10	100	60 18 121
U	2	PP	■	individually in bags 100 x 1	100	60 18 123
F	1	PS	–	bags 10 x 10	100	60 18 211
F	1	PS	■	individually in bags 100 x 1	100	60 18 213
F	2	PP	–	bags 10 x 10	100	60 18 221
F	2	PP	■	individually in bags 100 x 1	100	60 18 223
V	1	PS	–	bags 10 x 10	100	60 18 311
V	1	PS	■	individually in bags 100 x 1	100	60 18 313
V	2	PP	–	bags 10 x 10	100	60 18 321
V	2	PP	■	individually in bags 100 x 1	100	60 18 323

### Volume approx. 400 µl

Bottom	Fig. No.	Material	Sterilized	Color	Packaging	Quantity per Pack	Order No.
U	3	PP	–	white	bags 10 x 10	100	60 18 810
U	2	PP	–	transparent	bags 10 x 10	100	60 18 812
U	2	PP	■	transparent	individually in bags 100 x 1	100	60 18 813
U	3	PS	–	white	bags 10 x 10	100	60 18 820
U	1	PS	–	transparent	bags 10 x 10	100	60 18 822
U	1	PS	■	transparent	individually in bags 100 x 1	100	60 18 823

# Micro Test Plates, Spot Test Plates



## 96-Well Micro Test Plates, volume approx. 300 µl

Polypropylene (PP) or polystyrene (PS).

- with U-bottom or F-bottom
- SBS standard micro test plate format
- 8 x 12 wells alphanumerically coded
- wells with raised rims for reduced cross-contamination and effective sealing with Ratiolab cover films (see on page 69)
- precise design, perfectly shaped for use with automatic systems
- optimum space for bar code labels on the side walls
- space-saving and stackable
- surface untreated

Bottom	Material	Color	Packaging	Quantity per Pack	Order No.
U	PP	black	bags 10 x 10	100	60 18 710
U	PS	black	bags 10 x 10	100	60 18 711
U	PP	white	bags 10 x 10	100	60 18 712
U	PS	white	bags 10 x 10	100	60 18 713
F	PP	black	bags 10 x 10	100	60 18 717
F	PS	black	bags 10 x 10	100	60 18 718
F	PP	white	bags 10 x 10	100	60 18 719
F	PS	white	bags 10 x 10	100	60 18 720

## Lids for Micro Test Plates

Optional in highly transparent polystyrene (PS) or polypropylene (PP).

Material	Sterilized	Packaging	Quantity per Pack	Order No.
PS	–	bags 10 x 10	100	60 18 402
PS	■	individually in bags 100 x 1	100	60 18 401
PP	–	bags 10 x 10	100	60 18 422
PP	■	individually in bags 100 x 1	100	60 18 421

## Spot Test Plates

Polystyrene (PS), white.

- 3 x 4 wells
- L x W x H: 135 x 109 x 10 mm

Ø Mould mm	Packaging	Quantity per Pack	Order No.
20	bag 1 x 100	100	60 18 060