

Multi Purpose Beakers and Multi Purpose Bottles

Multi Purpose beakers with screw-caps

Polypropylene (PP).

For sampling, transport and storage of solid, pasty and liquid samples.

- with red screw-cap
- supplied in non-sterile or sterile, individually packed



Volume ml	Fig. No.	Material	Graduation and Labeling area	Sterilized	Height with cap, mm	Ø cap, outer, mm	Packaging	Quantity per Pack	Order No.
30	1	PP	—	—	47	43	cardboard box 1 x 1000	1000	59 21 173
30	1	PP	—	■	47	43	bags 500 x 1	500	59 21 174
30	2	PP	■	—	82	34	cardboard box 1 x 700	700	59 21 177
30	2	PP	■	■	82	34	bags 500 x 1	500	59 21 178
60	3	PP	■	—	66	43	cardboard box 1 x 600	600	59 21 175
60	3	PP	■	■	66	43	bags 500 x 1	500	59 21 176
120	4	PP	■	—	73	63	cardboard box 1 x 500	500	59 21 182
120	4	PP	■	■	73	63	bags 250 x 1	250	59 21 183
200	5	PP	—	—	80	64	cardboard box 1 x 200	200	59 21 185
200	5	PP	—	■	80	64	bags 150 x 1	150	59 21 184
200	6	PS	—	■	82	64	bags 150 x 1	150	59 21 188



Multi Purpose bottles

Polyethylene, blue translucent (PE).

- inclusive lid with nonwoven fabric liner
- with graduations, light protection
- wide neck diameter 80 mm
- Volume: 2.5 l

Product	Packaging	Quantity per Pack	Order No.
Multi purpose bottle 2.5 l with self-adhesive labels	cardboard box 1 x 60	60	59 22 250