


Multistand Plus Rack and Rack for Microcentrifuge Tubes



Multistand Plus Rack for sample tubes

Polypropylene (PP).

The multistand plus is suitable for sample tubes and you have the opportunity to store and place 50 tubes. The multistand plus has a temperature resistance up to 75°C.

- for tubes Ø 12 to Ø 18 mm
- autoclavable

Color	Positions	W x D x H mm	Quantity per Pack	Order No.
red	5 x 10	115 x 200 x 50	1	37 12 301
white	5 x 10	115 x 200 x 50	1	37 12 302

Rack for Microcentrifuge Tubes

Polypropylene (PP).

Ideal for the secure handling of reaction vessels. Space-saving on the laboratory bench due to its narrow form. Practical side grips.

- for tubes up to Ø 11 mm
- autoclavable at 121°C
- temperature-resistant from -80°C to +121°C
- 225 mm x 65 mm x 28 mm (LxWxH)



Color	For Micro Tubes	Positions	Quantity per Pack	Order No.
green	1.5 – 2.0 ml	5 x 16	5	37 12 006
pink	1.5 – 2.0 ml	5 x 16	5	37 12 008
blue	1.5 – 2.0 ml	5 x 16	5	37 12 005
orange	1.5 – 2.0 ml	5 x 16	5	37 12 007
yellow	1.5 – 2.0 ml	5 x 16	5	37 12 009
mix pack	1,5 – 2,0 ml	5 x 16	5	37 12 010