

Waste Disposal Bags, Waste Disposal Bags BIOHAZARD, Autoclavable at 121 °C

Custom sizes on request,
Autoclavable at 121 °C



Waste Disposal Bags, Dispenser Box, High Transparency

Polypropylene (PP).

Bags of very strong PP-film with a broad seam at the bottom, for the disposal of contaminated labware.

- autoclavable at 121 °C, 20 min, for autoclaving, the bag must not be completely closed
- for bag size 200 x 300 mm: practical bench-top stand, plastic-coated, with 3 rubber feet

| W x H mm | Volume approx. l | Film Thickness | Packaging | Quantity per Pack | Order No. |
|--|------------------|----------------|-------------------|-------------------|-----------|
| 200 x 300 | 1.5 | 50 µm | dispenser 1 x 100 | 100 | 70 01 005 |
| bench-top stand for waste disposal bags 200 x 300 mm, Height 25 cm | | | | 1 | 70 01 400 |
| Set consisting of bench-top stand and 100 waste disposal bags (high transparency) 200 x 300 mm | | | | 1 | 70 01 405 |



Waste Disposal Bags, High Transparency

Polypropylene (PP).

Bags of very strong PP-film with a broad seam at the bottom, for the disposal of contaminated labware.

- autoclavable at 121 °C, 20 min, for autoclaving, the bag must not be completely closed

| W x H mm | Volume approx. l | Film Thickness | Packaging | Quantity per Pack | Order No. |
|------------|------------------|----------------|-----------------------|-------------------|-----------|
| 300 x 500 | 6 | 50 µm | cardboard box 1 x 500 | 500 | 70 01 105 |
| 400 x 780 | 22 | 50 µm | cardboard box 1 x 500 | 500 | 70 01 155 |
| 600 x 800 | 32 | 50 µm | cardboard box 1 x 500 | 500 | 70 01 205 |
| 700 x 1100 | 72 | 50 µm | cardboard box 1 x 350 | 350 | 70 01 305 |



Waste Disposal Bags BIOHAZARD red – Warning Instruction and sterilization indicator patch

Polyethylene (PE).

Distinctive red autoclavable bags for easy identification, for the disposal of contaminated labware.

- rectangle white indicator patch field turns into dark when autoclaved
- with BIOHAZARD symbol and safety imprint in English, German, French and Spanish languages
- autoclavable at 121 °C, 20 min, for autoclaving, the bag must not be completely closed
- practical bench-top stand, plastic-coated, with 3 rubber feet
- only in large bag autoclaving

| W x H mm | Volume approx. l | Film Thickness | Packaging | Quantity per Pack | Order No. |
|--|------------------|----------------|------------------------|-------------------|-----------|
| 210 x 290 | 1.5 | 20 µm | cardboard box 10 x 100 | 1000 | 70 02 110 |
| Set consisting of bench-top stand and 100 red waste disposal bags 210 x 290 mm | | | | 1 | 70 01 410 |

Waste Disposal Bags BIOHAZARD, Autoclavable at 121 °C

Custom sizes on request,
Autoclavable at 121 °C



Waste Disposal Bags BIOHAZARD yellow – Warning Instruction and sterilization indicator patch

Polyethylene (HD-PE).

Autoclavable disposable waste bag of very strong PE-film for the sterilization and disposal of hazardous and contaminated labware.

- rectangle white indicator patch field turns into dark when autoclaved
- with BIOHAZARD symbol and safety imprint in German, English and French language
- autoclavable at 121 °C, 20 min, for autoclaving, the bag must not be completely closed

| W x H mm | Volume approx. l | Film Thickness | Packaging | Quantity per Pack | Order No. |
|------------|------------------|----------------|-----------------------|-------------------|-----------|
| 300 x 500 | 6 | 50 µm | cardboard box 1 x 500 | 500 | 70 03 100 |
| 400 x 780 | 22 | 50 µm | cardboard box 1 x 500 | 500 | 70 03 110 |
| 600 x 800 | 32 | 50 µm | cardboard box 1 x 500 | 500 | 70 03 120 |
| 700 x 1100 | 72 | 50 µm | cardboard box 1 x 350 | 350 | 70 03 130 |



Waste Disposal Bags, high transparency, BIOHAZARD

Polypropylene (PP).

For the sterilization and disposal of hazardous and contaminated labware.

- with blue BIOHAZARD symbol
- autoclavable at 121 °C, 20 min, for autoclaving, the bag must not be completely closed

| W x H mm | Volume | Film Thickness | Packaging | Quantity per Pack | Order No. |
|-----------|--------|----------------|----------------------|-------------------|-----------|
| 300 x 650 | 9 | 40 µm | cardboard box 4 x 50 | 200 | 70 02 111 |
| 415 x 580 | 30 | 40 µm | cardboard box 4 x 50 | 200 | 70 02 112 |
| 640 x 800 | 60 | 40 µm | cardboard box 4 x 50 | 200 | 70 02 113 |

Waste Disposal Bags, Waste Disposal Bags BIOHAZARD, Autoclavable at 134 °C

Custom sizes on request,
Autoclavable at 134 °C



Waste Disposal Bags, Dispenser Box, Standard

Polypropylene (PP).

Bags of very strong PP-film with a broad seam at the bottom, for the disposal of contaminated labware.

- autoclavable at 134 °C, 20 min, for autoclaving, the bag must not be completely closed
- for bag size 200 x 300 mm: practical bench-top stand, plastic-coated, with 3 rubber feet

| W x H mm | Volume approx. l | Film Thickness | Packaging | Quantity per Pack | Order No. |
|---|------------------|----------------|-------------------|-------------------|-----------|
| 200 x 300 | 1.5 | 50 µm | dispenser 1 x 100 | 100 | 70 01 000 |
| bench-top stand for waste disposal bags 200 x 300 mm, Height 25 cm | | | | 1 | 70 01 400 |
| Set consisting of bench-top stand and 100 waste disposal bags (standard) 200 x 300 mm | | | | 1 | 70 01 401 |

Waste Disposal Bags, Standard

Polypropylene (PP).

Bags of very strong PP-film with a broad seam at the bottom, for the disposal of contaminated labware.

- autoclavable at 134 °C, 20 min, for autoclaving, the bag must not be completely closed



| W x H mm | Volume approx. l | Film Thickness | Packaging | Quantity per Pack | Order No. |
|------------|------------------|----------------|-----------------------|-------------------|-----------|
| 300 x 500 | 6 | 50 µm | cardboard box 1 x 500 | 500 | 70 01 100 |
| 400 x 780 | 22 | 50 µm | cardboard box 1 x 500 | 500 | 70 01 150 |
| 600 x 800 | 32 | 50 µm | cardboard box 1 x 500 | 500 | 70 01 200 |
| 600 x 800 | 32 | 50 µm | cardboard box 1 x 100 | 100 | 70 01 210 |
| 600 x 800 | 32 | 100 µm | cardboard box 1 x 125 | 125 | 70 01 212 |
| 700 x 1100 | 72 | 50 µm | cardboard box 1 x 350 | 350 | 70 01 300 |
| 700 x 1100 | 72 | 50 µm | cardboard box 1 x 75 | 75 | 70 01 310 |
| 700 x 1100 | 72 | 100 µm | cardboard box 1 x 75 | 75 | 70 01 311 |

Waste Disposal bags BIOHAZARD – Warning Instruction and sterilization indicator patch

Polypropylene (PP).

Autoclavable disposable waste bag of very strong PP-film for the sterilization and disposal of hazardous and contaminated labware.

- rectangle white indicator patch field turns into dark when autoclaved
- with BIOHAZARD symbol and safety imprint in English, German, French and Spanish language
- autoclavable at 134 °C, 20 min, for autoclaving, the bag must not be completely closed



| W x H mm | Volume approx. l | Film Thickness | Packaging | Quantity per Pack | Order No. |
|------------|------------------|----------------|-----------------------|-------------------|-----------|
| 300 x 515 | 6 | 50 µm | cardboard box 1 x 100 | 100 | 70 02 105 |
| 300 x 515 | 6 | 50 µm | cardboard box 1 x 500 | 500 | 70 02 100 |
| 400 x 785 | 22 | 50 µm | cardboard box 1 x 100 | 100 | 70 02 155 |
| 400 x 785 | 22 | 50 µm | cardboard box 1 x 500 | 500 | 70 02 150 |
| 570 x 785 | 32 | 50 µm | cardboard box 1 x 100 | 100 | 70 02 205 |
| 570 x 785 | 32 | 50 µm | cardboard box 1 x 500 | 500 | 70 02 200 |
| 670 x 1085 | 72 | 50 µm | cardboard box 1 x 75 | 75 | 70 02 305 |
| 670 x 1085 | 72 | 50 µm | cardboard box 1 x 350 | 350 | 70 02 300 |

Waste Disposal Bags BIOHAZARD, Autoclavable at 134°C

Custom sizes on request,
Autoclavable at 134°C



Waste Disposal Bags BIOHAZARD red – Warning Instruction and sterilization indicator patch

Polypropylene (PP).

Distinctive red autoclavable bags for easy identification, for the disposal of contaminated labware.

- rectangle white indicator patch field turns into dark when autoclaved
- with BIOHAZARD symbol and safety imprint in German, English and French language
- autoclavable at 134°C, 20 min, for autoclaving, the bag must not be completely closed

| W x H mm | Volume approx. l | Film Thickness | Packaging | Quantity per Pack | Order No. |
|------------|------------------|----------------|-----------------------|-------------------|----------------------|
| 300 x 500 | 6 | 50 µm | cardboard box 1 x 500 | 500 | 70 03 200 |
| 400 x 780 | 22 | 50 µm | cardboard box 1 x 500 | 500 | 70 03 210 |
| 600 x 800 | 32 | 50 µm | cardboard box 1 x 500 | 500 | 70 03 220 |
| 600 x 800 | 32 | 50 µm | cardboard box 1 x 100 | 100 | New 70 03 225 |
| 700 x 1100 | 72 | 50 µm | cardboard box 1 x 350 | 350 | 70 03 230 |
| 700 x 1100 | 72 | 50 µm | cardboard box 1 x 75 | 75 | New 70 03 235 |

Waste Disposal Bags BIOHAZARD yellow – Warning Instruction and sterilization indicator patch

Polypropylene (PP).

Autoclavable disposable waste bag of very strong PP-film for the sterilization and disposal of hazardous and contaminated labware.

- rectangle white indicator patch field turns into dark when autoclaved
- with BIOHAZARD symbol and safety imprint in German, English and French language
- autoclavable at 134°C, 20 min, for autoclaving, the bag must not be completely closed



| W x H mm | Volume approx. l | Film Thickness | Packaging | Quantity per Pack | Order No. |
|------------|------------------|----------------|-----------------------|-------------------|----------------------|
| 300 x 500 | 6 | 50 µm | cardboard box 1 x 500 | 500 | New 70 03 300 |
| 400 x 780 | 22 | 50 µm | cardboard box 1 x 500 | 500 | New 70 03 310 |
| 600 x 800 | 32 | 50 µm | cardboard box 1 x 500 | 500 | New 70 03 320 |
| 600 x 800 | 32 | 50 µm | cardboard box 1 x 100 | 100 | New 70 03 325 |
| 700 x 1100 | 72 | 50 µm | cardboard box 1 x 350 | 350 | New 70 03 330 |
| 700 x 1100 | 72 | 50 µm | cardboard box 1 x 75 | 75 | New 70 03 335 |

Bench-top Stand, Transport bag BIOHAZARD and Cleanibag System

Bench-top stand for waste disposal bags

Practical bench-top stand, plastic-coated, with 3 rubber feet.

- for bag size 200 x 300 mm and 300 x 500 mm



| Product | Bild-Nr. | Quantity per Pack | Order No. |
|--|----------|-------------------|-------------------------|
| Bench-top stand for waste disposal bags 200 x 300 mm, Height 25 cm | 2 | 1 | 70 01 400 |
| Set consisting of bench-top stand and 100 waste disposal bags (standard) 200 x 300 mm | 3 | 1 | 70 01 401 |
| Set consisting of bench-top stand and 100 waste disposal bags (high transparency) 200 x 300 mm | 3 | 1 | 70 01 405 |
| bench-top stand for waste disposal bags 300 x 500 mm, Height 35 cm | 1 | 1 | 70 01 430 |
| Set consisting of bench-top stand and 100 waste disposal bags (standard) 300 x 500 mm | 3 | 1 | 70 01 415 |
| Set consisting of bench-top stand and 100 waste disposal bags (high transparency) 300 x 500 mm | 3 | 1 | 70 01 420 |
| Set consisting of bench-top stand and 100 red waste disposal bags 210 x 290 mm | 4 | 1 | 70 01 410 |
| Bench top stand with tilt base, for waste disposal bags 200 x 300 mm | 5 | 1 | New 70 01 440 |
| Design table stand made of stainless steel autoclavable, for waste disposal bags 200 x 300 mm | 6 | 1 | New 70 01 450 |



Transport bags BIOHAZARD

Polyethylene (LD-PE).

- transparent with fast-zip closure for reclosing and with a separate external pouch for documentation
- suitable for sample transport

| W x H mm | Film Thickness | Packaging | Quantity per Pack | Order No. |
|-----------------|----------------|-----------------------|-------------------|-----------|
| 155 x 250 / 210 | 50 µm | cardboard box 1 x 100 | 100 | 70 00 300 |
| 230 x 305 / 250 | 50 µm | cardboard box 1 x 100 | 100 | 70 00 310 |



Cleanibag System

All-purpose PE bags with self-adhesive strip for clinical and laboratory use.

- self-adhesive strip, sticks on smooth and clean surfaces
- residue-free on removable

Other sizes available on request.

| W x H mm | Volume approx. l | Packaging | Quantity per Pack | Order No. |
|-----------|------------------|------------------------|-------------------|-----------|
| 200 x 300 | 1.5 | cardboard box 1 x 1000 | 1000 | 70 00 250 |
| 300 x 300 | 6.5 | cardboard box 1 x 1000 | 1000 | 70 00 260 |
| 300 x 450 | 8 | cardboard box 1 x 1000 | 1000 | 70 00 270 |

Self-adhesive