

Cryo Boxes, Cardboard 1/2-Format and Cryo Labels



Cryo Boxes and Grid Inserts, Cardboard, 1/2-Format 130 x 62 mm

Ultra-low temperature-resistant cardboard, white.

- for sample storage even at ultra-low temperatures
- temperature-resistant down to -80°C
- water-repellent finish

Product	Grid	W x D x H mm	For Tube Dia. mm	Quantity per Pack	Order No.
1/2-format Cryo Boxes, cardboard	–	130 x 62 x 50	–	10	51 31 506
1/2-format Grid Inserts	8 x 4	124 x 58 x 30	14.0	10	51 20 218
1/2-format Grid Inserts	10 x 5	124 x 58 x 30	10.5	10	51 20 128



Rainbow Cryo labels

For customizing the white Cryo Boxes, the Rainbow Cryo labels are available in modern colors: red, yellow, green, blue, purple.

- temperature-resistant down to -150°C
- 50 labels per pack, divided into 5 colors, 10 pieces from each color

Product	L x W mm	Quantity per Pack	Order No.
Rainbow Cryo labels	200 x 40	5 x 10	51 20 002