

# Micro Centrifuge and Thermocult Incubator



## Micro Centrifuge for PCR Strips and Microcentrifuge Tubes

Space- and energy-saving micro centrifuge with adequate power for routine centrifugations.

- speed 7000/min, max RCF 2680 x g
- lid must be closed for operation mode
- safety: the lid can only be opened when the rotor is not rotating
- rapid start and stop within seconds
- power supply AC100-240V/50Hz/60Hz
- diameter 150 mm, height 117 mm, weight 0.5 kg
- 1 year guarantee
- supplied complete with two rotors:
  - rotor 1 for: 8 microcentrifuge tubes 0.5 ml till 2.0ml
  - rotor 2 for: 16 microcentrifuge tubes 0.2 ml or two pcr strips of 8 tubes 0.2 ml

Product	Quantity per Pack	Order No.
Micro Centrifuge with two rotors	1	71 00 007

## Thermocult Incubator

Compact incubator for bacteriology. Fixed temperature setting at  $35 \pm 2^\circ\text{C}$ , adjustable up to  $42^\circ\text{C}$ , with built-in control thermometer. Glass-clear plastic door for visual control. The incubator can be placed on work benches or set up on walls.



- inside size (W x D x H): 210 x 145 x 110 mm, for up to 20 tissue cultures or 8 petri dishes
- outside size (W x D x H): 300 x 170 x 150 mm
- weight: approx. 1.3 kg, power: 230 V, 30 W
- shelf available separately
- 2 year guarantee

Product	Quantity per Pack	Order No.
Thermocult Incubator	1	90 22 000
shelf	1	90 22 001