


# Cryo Boxes with Slip Lid, PP, Format 133 x 133 mm



## 2-in-1 Cryo Boxes with Slip Lid, PP, Format 133 x 133 mm


High-quality polypropylene (PP). Available in natural and other colors.

The specialty: the boxes are adjustable in height.

The boxes are available in 6 modern colors, with numerical coding and corresponding vents for draining condensed water.

The 6 x 6 until 10 x 10 grid accommodates space from 36 to 100 tubes.

- autoclavable, temperature-resistant from  $-90^{\circ}\text{C}$  to  $+121^{\circ}\text{C}$
- grids are available: 6 x 6, 7 x 7, 8 x 8, 9 x 9, 10 x 10
- numerically coded



	Grid For tube Dia.	6 x 6 19.0 mm	7 x 7 17.0 mm	8 x 8 15.0 mm	9 x 9 13.0 mm	10 x 10 11.5 mm
Color	Quantity per Pack	Order No.	Order No.	Order No.	Order No.	Order No.
natural	5	53 20 410	53 20 310	53 20 210	53 20 010	53 20 110
yellow	5	53 20 411	53 20 311	53 20 211	53 20 011	53 20 111
red	5	53 20 412	53 20 312	53 20 212	53 20 012	53 20 112
green	5	53 20 413	53 20 313	53 20 213	53 20 013	53 20 113
blue	5	53 20 414	53 20 314	53 20 214	53 20 014	53 20 114
black	5	53 20 415	53 20 315	53 20 215	53 20 015	53 20 115
mixed	5	53 20 416	53 20 316	53 20 216	53 20 016	53 20 116

# Cryo Boxes with Grids, PP, Format 133 x 133 mm




## Cryo Boxes with Grids, PP, Format 133 x 133 mm

The boxes are made of high quality polypropylene in natural or other colors. They are available with heights of 52 or 75 mm.

- autoclavable, temperature-resistant from – 90 °C to + 121 °C
- 81 places, 9 x 9 grid, numerically coded
- robust hinge with safe snap-on lid
- stackable for safe transport
- drain bores at the bottom for dew liquids
- removable grid inserts
- Cryo Boxes can also be used for larger vials by removing the grid inserts

9 x 9 grid,  
numerically coded



Color	W x D x H mm	For Tube Dia. mm	Grid	Quantity per Pack	Order No.
<b>for Cryo Tubes 1.2 – 2 ml</b>					
natural	133 x 133 x 52	12.8	9 x 9	5	51 20 010
yellow	133 x 133 x 52	12.8	9 x 9	5	51 20 011
red	133 x 133 x 52	12.8	9 x 9	5	51 20 012
green	133 x 133 x 52	12.8	9 x 9	5	51 20 013
blue	133 x 133 x 52	12.8	9 x 9	5	51 20 014
black	133 x 133 x 52	12.8	9 x 9	5	51 20 015
<b>for Cryo Tubes 3 – 4 ml</b>					
natural	133 x 133 x 75	12.8	9 x 9	5	51 20 020
yellow	133 x 133 x 75	12.8	9 x 9	5	51 20 021
red	133 x 133 x 75	12.8	9 x 9	5	51 20 022
green	133 x 133 x 75	12.8	9 x 9	5	51 20 023
blue	133 x 133 x 75	12.8	9 x 9	5	51 20 024
black	133 x 133 x 75	12.8	9 x 9	5	51 20 025

# Cryo Boxes and Grid Inserts, PP, Format 133 x 133 mm



## Cryo Boxes and Grid Inserts, PP, Format 133 x 133 mm

The boxes are made of high quality polypropylene in natural or other colors. They are available with heights of 52 or 75 mm. These ratiolab® Cryo Boxes are particularly practical and flexible due to their removable grid inserts.

The separate naturally-colored inserts are available in five sizes and can thus be used with sample vessels of differing diameters. The inserts fit all ratiolab® Cryo Boxes. Cleaning of the boxes and the inserts is much simpler than with Cryo Boxes with permanently fitted inserts.

- autoclavable, temperature-resistant from – 90 °C to + 121 °C
- robust hinge with safe snap-on lid
- stackable for safe transport
- drain bores at the bottom for dew liquids
- Cryo Boxes can also be used for larger vials by removing the grid inserts

Color	W x D x H mm	For Tube Dia. mm	Quantity per Pack	Order No.
<b>Cryo Boxes without Grid Insert</b>				
natural	133 x 133 x 52		5	52 20 010
yellow	133 x 133 x 52		5	52 20 011
red	133 x 133 x 52		5	52 20 012
green	133 x 133 x 52		5	52 20 013
blue	133 x 133 x 52		5	52 20 014
black	133 x 133 x 52		5	52 20 015
natural	133 x 133 x 75		5	52 20 020
yellow	133 x 133 x 75		5	52 20 021
red	133 x 133 x 75		5	52 20 022
green	133 x 133 x 75		5	52 20 023
blue	133 x 133 x 75		5	52 20 024
black	133 x 133 x 75		5	52 20 025
<b>Grid Inserts</b>				
6 x 6	133 x 133 x 30	20.5	5	52 20 100
7 x 7	133 x 133 x 30	17.3	5	52 20 101
8 x 8	133 x 133 x 30	15.8	5	52 20 102
9 x 9	133 x 133 x 30	12.8	5	52 20 103
10 x 10	133 x 133 x 30	11.2	5	52 20 104