


Cryo Boxes with Grids, PC, Format 132 x 132 mm and 76 x 76 mm



Cryo Boxes with Grids, PC, Format 132 x 132 mm and 76 x 76 mm

Made of high quality colored polycarbonate (PC).
For use with liquid nitrogen.

- autoclavable, temperature-resistant from – 190 °C to + 121 °C
- 81 or 100 places, 9 x 9 resp. 10 x 10 grids
- 1/4-format with 25 places, 5 x 5 grid
- numerically coded
- stackable for safe transport
- drain bores at the bottom for dew liquids
- Cryo Boxes can also be used for larger vials by removing the grid inserts (except figure 2)

Color	Fig No.	W x D x H mm	For Tube Dia. mm	Grid	Quantity per Pack	Order No.
for Cryo Tubes 1.2 – 2 ml						
blue	1	132 x 132 x 52	12,3	9 x 9	5	51 50 001
green	1	132 x 132 x 52	12,3	9 x 9	5	51 50 002
red	1	132 x 132 x 52	12,3	9 x 9	5	51 50 003
yellow	1	132 x 132 x 52	12,3	9 x 9	5	51 50 004
for Cryo Tubes 1.2 – 2 ml						
blue	2	132 x 132 x 52	12,3	10 x 10	5	51 50 011
green	2	132 x 132 x 52	12,3	10 x 10	5	51 50 012
red	2	132 x 132 x 52	12,3	10 x 10	5	51 50 013
yellow	2	132 x 132 x 52	12,3	10 x 10	5	51 50 014
white	2	132 x 132 x 52	12,3	10 x 10	5	51 50 015
for Cryo Tubes 3 – 5 ml						
blue	3	132 x 132 x 94	12,3	9 x 9	4	51 50 021
green	3	132 x 132 x 94	12,3	9 x 9	4	51 50 022
red	3	132 x 132 x 94	12,3	9 x 9	4	51 50 023
yellow	3	132 x 132 x 94	12,3	9 x 9	4	51 50 024
1/4-format, for Cryo Tubes 1.2 – 2 ml						
blue	4	76 x 76 x 52	12,3	5 x 5	8	51 50 031
green	4	76 x 76 x 52	12,3	5 x 5	8	51 50 032
red	4	76 x 76 x 52	12,3	5 x 5	8	51 50 033
yellow	4	76 x 76 x 52	12,3	5 x 5	8	51 50 034