

# Racks for Freezer Chests



## Vertical Racks, stainless steel

For the storing of Cryo Boxes in freezer chests.  
Stainless steel.


- for all Cryo Boxes format 130 x 130; 133 x 133 and 136 x 136 mm
- surfaces for perfect cleaning
- robust construction
- handles for convenient transport
- rack size includes the holder
- Cryo Racks can be supplied with shelf sizes for any standard Cryo Boxes
- locking rod on request

Rack standard sizes are listed in the tables below:

other sizes  
available on request

Number of Shelves horizontal x vertical	Rack Size D x W x H mm	for Cryo Boxes Height			Quantity per Pack	Order No.
		32 mm	50 mm	75 mm		
1 x 4	140 x 140 x 222	–	■	–	1	54 00 041
1 x 5	140 x 140 x 277	–	■	–	1	54 00 021
1 x 6	140 x 140 x 482	–	–	■	1	54 00 028
1 x 7	140 x 140 x 385	–	■	–	1	54 00 042
1 x 7	140 x 140 x 561	–	–	■	1	54 00 029
1 x 8	140 x 140 x 441	–	■	–	1	54 00 027
1 x 8	140 x 140 x 641	–	–	■	1	54 00 030
1 x 9	140 x 140 x 495	–	■	–	1	54 00 026
1 x 9	140 x 140 x 720	–	–	■	1	54 00 031
1 x 10	140 x 140 x 550	–	■	–	1	54 00 000
1 x 11	140 x 140 x 406	■	–	–	1	54 00 032
1 x 11	140 x 140 x 604	–	■	–	1	54 00 023
1 x 12	140 x 140 x 443	■	–	–	1	54 00 033
1 x 12	140 x 140 x 659	–	■	–	1	54 00 025
1 x 13	140 x 140 x 479	■	–	–	1	54 00 034
1 x 13	140 x 140 x 713	–	■	–	1	54 00 001
1 x 14	140 x 140 x 516	■	–	–	1	54 00 035
1 x 15	140 x 140 x 552	■	–	–	1	54 00 036
1 x 16	140 x 140 x 589	■	–	–	1	54 00 037
1 x 17	140 x 140 x 625	■	–	–	1	54 00 038
1 x 18	140 x 140 x 662	■	–	–	1	54 00 024

# Racks for Freezer Cabinets



## Horizontal Racks, stainless steel

For the storing of Cryo Boxes in freezer cabinets and chests.  
Stainless steel.

- for all Cryo Boxes format 130 x 130; 133 x 133 and 136 x 136 mm
- surfaces for perfect cleaning
- robust construction
- handles for convenient transport
- Cryo Racks can be supplied with shelf sizes for any standard Cryo Boxes

Rack standard sizes are listed in the tables below:

**other sizes  
available on request**

Number of Shelves horizontal x vertical	Number of Boxes	Rack Size D x W x H mm	for Cryo Boxes Height				Quantity per Pack	Order No.
			50 mm	75 mm	100 mm	130 mm		
3 x 3	9	424 x 140 x 165	■	–	–	–	1	54 00 007
3 x 4	12	424 x 140 x 219	■	–	–	–	1	54 00 009
3 x 5	15	424 x 140 x 274	■	–	–	–	1	54 00 022
4 x 2	8	566 x 140 x 164	–	■	–	–	1	54 00 011
4 x 2	8	566 x 140 x 219	–	–	■	–	1	54 00 010
4 x 2	8	566 x 140 x 275	–	–	–	■	1	54 00 012
4 x 3	12	566 x 140 x 240	–	■	–	–	1	54 00 006
4 x 3	12	566 x 140 x 314	–	–	■	–	1	54 00 040
4 x 4	16	566 x 140 x 219	■	–	–	–	1	54 00 004
4 x 4	16	566 x 140 x 319	–	■	–	–	1	54 00 005
4 x 5	20	566 x 140 x 274	■	–	–	–	1	54 00 003
4 x 6	24	566 x 140 x 328	■	–	–	–	1	54 00 002
4 x 7	28	566 x 140 x 391	■	–	–	–	1	54 00 008
5 x 2	10	705 x 140 x 161	–	■	–	–	1	54 00 016
5 x 2	10	705 x 140 x 212	–	–	■	–	1	54 00 020
5 x 3	15	705 x 140 x 240	–	■	–	–	1	54 00 017
5 x 4	20	705 x 140 x 219	■	–	–	–	1	54 00 013
5 x 4	20	705 x 140 x 319	–	■	–	–	1	54 00 018
5 x 5	25	705 x 140 x 274	■	–	–	–	1	54 00 014
5 x 6	30	705 x 140 x 328	■	–	–	–	1	54 00 015

# Racks for Freezer Cabinets



## Drawer Rack, stainless steel

For the storing of Cryo Boxes in freezer cabinets and chests.  
Stainless steel.

- for all Cryo Boxes format 130 x 130; 133 x 133 and 136 x 136 mm
- cabinet slide-ins without pull-out stop
- easy-to-clean polished surface
- robust design

Rack standard sizes are listed in the tables below:

**other sizes  
available on request**

Number of Shelves horizontal x vertical	Number of Shelves	Rack Size D x W x H mm	for Cryo Boxes Height		Quantity per Pack	Order No.
			50 mm	75 mm		
4 x 2	8	570 x 144 x 112	■	–	1	54 00 043
4 x 3	12	570 x 144 x 168	■	–	1	54 00 044
4 x 3	12	570 x 144 x 250	–	■	1	54 00 019
4 x 4	16	570 x 144 x 224	■	–	1	54 00 046
4 x 4	16	570 x 144 x 335	–	■	1	54 00 066
4 x 5	20	570 x 144 x 280	■	–	1	54 00 045
4 x 6	24	570 x 144 x 336	■	–	1	54 00 047
4 x 7	28	570 x 144 x 392	■	–	1	54 00 049
5 x 3	15	709 x 144 x 241	–	■	1	54 00 052
5 x 4	20	709 x 144 x 224	■	–	1	54 00 050
5 x 4	20	709 x 144 x 321	–	■	1	54 00 053
5 x 5	25	709 x 144 x 280	■	–	1	54 00 048
5 x 6	30	709 x 144 x 336	■	–	1	54 00 051