

free of RNase-, DNase-, DNA, ATP-,
pyrogens, and PCR inhibitors-,
Low Retention

ratiolab[®]
disposables for sciences

Product information
lot number and expiry date

Practical hinged lid
opens up to 180°

Pipet tips and aerosol filters
made of ultra-hydrophobic HD-PE

Safe click-closure
for dust protection

Stackable and space-saving



Liquid handling Products



THINK! Development and production
of customized products



Liquid handling
products



Labware for general
lab applications



Life science products for
microbiology and biotechnology



Cryo products for
deep-freeze applications



Slamed Pipets and
Analysers

aeroject[®] ultra filter tips – Pipet Tips with Aerosol Barrier

ratiolab[®] pipet tips are ideal for all commercially available pipets. Their uniform high quality and precision make today's pipets fully effective.

aeroject[®] tips with special filters made of hydrophobic HD-PE as an aerosol barrier and for protection against contamination.

Ideal fit for manual or electronic single- and multichannel pipettes.

Suitable for all commonly used pipets.





aeroject® ultra filter tips – Pipet Tips with Aerosol Barrier



Pipet tips with special filters, to effectively prevent the pipettor shaft from contamination with aerosols and liquids to be sucked into the pipettor.

- sterilized
- free of RNase, DNase, pyrogens and ATP
- free of human DNA and PCR inhibitors
- pipet tips and aerosol filters made of ultra-hydrophobic HD-PE
- minimised retention due to special surface treatment
- perfect and tight fit on all major pipettor brands
- aeroject® ultra 17 70 010 and 17 71 000xl with extended length for operation with narrow vessels and reduced risk of cross-contamination
- packed in plastic racks with hinged lids

free of RNase-, DNase-, DNA, ATP-,
pyrogens, and PCR inhibitors-,
Low Retention

Type	Volume Range μ l	Length mm	Packaging	Quantity per Pack	Order No.
Pipet Tip aeroject® ultra tips xl, for all commonly used pipets, ratiopetta®	0,1-10	45,60	Racks 10 x 96	960	17 70 010
					
Pipet Tip aeroject® ultra tips, for all commonly used pipets, ratiopetta®	1-20	50,80	Racks 10 x 96	960	17 70 020
					
Pipet Tip aeroject® ultra tips, for all commonly used pipets, ratiopetta®	1-100	50,80	Racks 10 x 96	960	17 70 100
					
Pipet Tip aeroject® ultra tips, for all commonly used pipets, ratiopetta®	1-200	50,80	Racks 10 x 96	960	17 70 200
					
Pipet Tip aeroject® ultra tips, for all commonly used pipets, ratiopetta®	100-1000	86,00	Racks 10 x 96	960	17 71 000xl
